

Product Datasheet

PAX3 Antibody (orb751478)

Description	PAX3 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Peptide with sequence TTLAGAVPRMM-C, from the N Terminus of the protein sequence according to NP_000429.2; NP_039230.1; NP_852122.1;.1 NP_852123.1; NP_852124.1; NP_852125.1; NP_852126.1.
Form/Appearance	Affinity Purified
Storage	Antibody can be aliquotted and stored frozen at -20°; C to -70°; C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2°; - 8°; C for 1 month without detectable loss of activity. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Goat IgG
Uniprot ID	P23760
Entrez	5077
Dilution Range	IHC 2-15 ug/mL; WB 0.1-0.2 ug/mL; E 0.5-1.0 ug/mL
Expiration Date	12 months from date of receipt.