



Product Datasheet

TMEM18 Antibody (orb75148)

Description	TMEM18 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96B42
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of rat brain lysat...