

Product Datasheet

Tau (Tau 46) Antibody (orb751485)

Description	Tau (Tau 46) Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Native bovine tau and the epitope maps to the carboxy-terminus of the protein.
Form/Appearance	Protein G Purified
Storage	Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P10636
Entrez	4137
Expiration Date	12 months from date of receipt.