

## Product Datasheet

### mCherry Antibody (orb751490)

|                            |                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | mCherry Antibody                                                                                            |
| <b>Species/Host</b>        | Mouse                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                                                                            |
| <b>Immunogen</b>           | Full length recombinant mCherry protein                                                                     |
| <b>Form/Appearance</b>     | Affinity Purified                                                                                           |
| <b>Storage</b>             | Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.                                |
| <b>Note</b>                | For research use only                                                                                       |
| <b>Isotype</b>             | Mouse IgG                                                                                                   |
| <b>Dilution Range</b>      | Immunohistochemistry: 1:500<br>Immunofluorescent: 1:500<br>Immunocytochemistry: 1:500 Western Blot: 1:2,000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                             |