

Product Datasheet

PEA-15 Antibody (orb751495)

Description	PEA-15 Antibody
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Full length human PEA-15 as expressed in and purified from E. coli
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	Q15121
Entrez	8682
Dilution Range	Immunofluorescent: 1:1,000-2,000 Immunocytochemistry: 1:1,000-2,000 Western Blot: 1:1,000-2,000
Expiration Date	12 months from date of receipt.