

Product Datasheet

FGF-2 Antibody (orb751499)

Description	FGF-2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Human FGF2
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgM
Uniprot ID	P09038
Entrez	2247
Expiration Date	12 months from date of receipt.