

Product Datasheet

PSAP (Prostatic Specific Acid Phosphatase)/PAP Antibody (orb751559)

Description	PSAP (Prostatic Specific Acid Phosphatase)/PAP Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant Human PSAP/Prosaposin
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P15309
Entrez	55
Dilution Range	Immunohistochemistry: 1:1,000-3,000
Expiration Date	12 months from date of receipt.