

Product Datasheet

Presenilin 2 Antibody (orb751613)

Description	Presenilin 2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Peptide corresponding to residues surrounding amino acid 330 of human presenilin 2
Form/Appearance	Affinity Purified
Storage	Do not aliquot. Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P49810
Entrez	5664
Expiration Date	12 months from date of receipt.