

Product Datasheet

phospho-NMDA Receptor 1 (Ser890) Antibody (orb751619)

Description	phospho-NMDA Receptor 1 (Ser890) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Phosphopeptide corresponding to residues around serine 890 of human NMDAR1
Form/Appearance	Affinity Purified
Storage	Do not aliquot. Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q05586
Entrez	2902
Expiration Date	12 months from date of receipt.