

Product Datasheet

Parkin Antibody (orb751621)

Description	Parkin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Synthetic peptide corresponding to residues around amino acid 400 of human parkin
Form/Appearance	Affinity Purified
Storage	Do not aliquot. Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	O60260
Entrez	5071
Expiration Date	12 months from date of receipt.