

Product Datasheet

c-Jun Antibody (orb751626)

Description	c-Jun Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Peptide derived from the amino-terminal sequence of human c-Jun
Form/Appearance	Affinity Purified
Storage	Do not aliquot. Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P05412
Entrez	3725
Expiration Date	12 months from date of receipt.