

Product Datasheet

alpha-Synuclein Antibody (orb751628)

Description	alpha-Synuclein Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Peptide derived from the carboxy-terminal sequence of human alpha-synuclein.
Form/Appearance	Affinity Purified
Storage	Do not aliquot. Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P37840
Entrez	6622
Expiration Date	12 months from date of receipt.