

Product Datasheet

phospho-Tyrosine Hydroxylase (Ser40) Antibody (orb751630)

Description	phospho-Tyrosine Hydroxylase (Ser40) Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Synthetic phospho-peptide (KLH coupled) derived from the sequence surrounding Ser40 of human tyrosine hydroxylase.
Form/Appearance	Affinity Purified
Storage	Do not aliquot. Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P07101
Entrez	7054
Expiration Date	12 months from date of receipt.