

Product Datasheet

SARS-CoV2 S-Protein ACE2 Binding Domain Antibody (orb751639)

Description	SARS-CoV2 S-Protein ACE2 Binding Domain Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant SARS-CoV2 S-Protein ACE2 binding domain expressed in and purified from E. coli
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P0DTC2
Dilution Range	Immunofluorescent: 1:3,000-5,000 Immunocytochemistry: 1:3,000-5,000 Western Blot: 1:3,000-5,000
Expiration Date	12 months from date of receipt.