



Product Datasheet

ATOH8 Antibody (orb75173)

Description	ATOH8 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96SQ7
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of A-20 cell lysat...