

Product Datasheet

IL-4R/CD124 (phospho-Y497) antibody (orb751923)

Description	Rabbit polyclonal antibody to IL-4R/CD124 (phospho-Y497)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human IL-4R/CD124 around the phosphorylation site of Tyrosine 497.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.4.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 78 kDa
Uniprot ID	P24394
Dilution Range	WB: 1: 500-1000
Expiration Date	12 months from date of receipt.