



Product Datasheet

erbB-2 Antibody (orb75194)

Description	erbB-2 Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04626
Expiration Date	12 months from date of receipt.



Western blot analysis of rat liver tissue...



Immunohistochemical analysis of formalin-fixed tissue...