

## Product Datasheet

### MCL1 (phospho-Ser159) antibody (orb751965)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MCL1 (phospho-Ser159)                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | ELISA, IF, WB                                                                                                                              |
| <b>Immunogen</b>           | The antiserum was produced against synthesized peptide derived from human MCL1 around the phosphorylation site of Ser159. AA range:125-174 |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.                                                               |
| <b>Concentration</b>       | 1mg/ml                                                                                                                                     |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.                                                   |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                 |
| <b>MW</b>                  | ~ 39 kDa                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>Q07820</b>                                                                                                                              |
| <b>Dilution Range</b>      | IF: 1: 200-1000, ELISA: 1: 10000                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |