

Product Datasheet

Trk A (phospho-Tyr496) antibody (orb751968)

Description	Rabbit polyclonal antibody to Trk A (phospho-Tyr496)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Trk A around the phosphorylation site of Tyr496. AA range:471-520
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	P04629
Dilution Range	IF: 1: 200-1000, ELISA: 1: 5000
Expiration Date	12 months from date of receipt.