



Product Datasheet

LDL-R Antibody (orb75199)

Description	LDL-R Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01130
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver tis...



Immunohistochemical analysis of formalin...