

Product Datasheet

NMDAR2B (phospho-Tyr1336) antibody (orb752006)

Description	Rabbit polyclonal antibody to NMDAR2B (phospho-Tyr1336)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human NMDAR2B around the phosphorylation site of Tyr1336. AA range:1302-1351
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 150 kDa
Uniprot ID	Q13224
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, ELISA: 1: 10000
Expiration Date	12 months from date of receipt.