

Product Datasheet

SH2D1B antibody (orb752010)

Description	Rabbit polyclonal antibody to SH2D1B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	The antiserum was produced against synthesized peptide derived from human SH2D1B. AA range:71-120
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	/
Uniprot ID	O14796
Dilution Range	ICC: 1: 100-300, ELISA: 1: 40000
Expiration Date	12 months from date of receipt.