

Product Datasheet

RECK antibody (orb752078)

Description	Rabbit polyclonal antibody to RECK
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Synthetic peptide, corresponding to Human RECK.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 110 kDa
Uniprot ID	O95980
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, IF: 1: 200-1000, ELISA: 1: 20000
Expiration Date	12 months from date of receipt.