

Product Datasheet

CA1 antibody (orb752114)

Description	Rabbit polyclonal antibody to CA1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CA1. AA range:41-90
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 29 kDa
Uniprot ID	P00915
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, ELISA: 1: 40000
Expiration Date	12 months from date of receipt.