

Product Datasheet

IFI16 antibody (orb752124)

Description	Rabbit polyclonal antibody to IFI16
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human IFI16. AA range:731-780
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 88 kDa
Uniprot ID	Q16666
Dilution Range	WB: 1: 500-2000, ELISA: 1: 40000
Expiration Date	12 months from date of receipt.