

Product Datasheet

PC-PLD2 (phospho Tyr169) antibody (orb752215)

Description	Rabbit polyclonal antibody to PC-PLD2 (phospho Tyr169)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PLD2 around the phosphorylation site of Tyr169. AA range:136-185
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~95 kDa
Uniprot ID	O14939
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, ELISA: 1: 10000
Expiration Date	12 months from date of receipt.