

Product Datasheet

CD227/MUC1 (phospho-Tyr1229) antibody (orb752239)

Description	Rabbit polyclonal antibody to CD227/MUC1 (phospho-Tyr1229)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CD227/MUC1 around the phosphorylation site of Tyr1229. AA range:1201-1250
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 170 kDa
Uniprot ID	P15941
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, IF: 1: 200-1000, ELISA: 1: 10000
Expiration Date	12 months from date of receipt.