

Product Datasheet

hnRPD (phospho-Ser83) antibody (orb752240)

Description	Rabbit polyclonal antibody to hnRPD (phospho-Ser83)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human hnRPD around the phosphorylation site of Ser83. AA range:49-98
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 38 kDa
Uniprot ID	Q14103
Dilution Range	WB: 1: 500-2000, ICC: 1: 100-300, IF: 1: 200-1000, ELISA: 1: 20000
Expiration Date	12 months from date of receipt.