

Product Datasheet

APOB(6G3-C7) monoclonal antibody (orb752258)

Description	APOB(6G3-C7) monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant full length Human APOB .
Preservatives	Mouse IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Storage	APOB monoclonal antibody detects endogenous levels of APOB protein.
Note	For research use only
Clonality	Monoclonal
MW	~ 515kDa
Uniprot ID	P04114
Dilution Range	ELISA WB
Expiration Date	12 months from date of receipt.