


Product Datasheet

Swine H1N1 Neuraminidase Antibody (orb75231)

Description	Swine H1N1 Neuraminidase Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	C4AL28
Expiration Date	12 months from date of receipt.



Bar chart illustrates about the ELISA re...