


Product Datasheet

Ang4 antibody (orb752360)

Description	Ang4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Angiogenin-4 protein (25-144aa)
Target	Ang4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3TMQ6
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.



WB analysis of Mmouse heart tissue, mouse...