


Product Datasheet

ompF antibody (orb752372)

Description	ompF antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli Outer membrane protein F protein (23-362aa)
Target	OMPf
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 95%, Protein G purified
Uniprot ID	P02931
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.



WB detected
in
Recombinant
protein
using...