

Product Datasheet

AGO3 antibody (orb752384)

Description	AGO3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Protein argonaute-3 protein (251-400aa)
Target	AGO3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H9G7
Dilution Range	IHC-P: 1:100-1:500
Expiration Date	12 months from date of receipt.