

Product Datasheet

AGO4 antibody (orb752388)

Description	AGO4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Protein argonaute-4 protein (716-861aa)
Target	AGO4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCK5
Dilution Range	IF: 1:100-1:500
Expiration Date	12 months from date of receipt.