

Product Datasheet

ATAD3A antibody (orb752400)

Description	ATAD3A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human ATPase family aaA domain-containing protein 3A protein (2-240aa)
Target	ATAD3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NVI7
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.