

Product Datasheet

C1QTNF9 antibody (orb752412)

Description	C1QTNF9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Complement C1q and tumor necrosis factor-related protein 9A protein (58-72aa)
Target	C1QTNF9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C862
Dilution Range	IHC-P: 1:100-1:500
Expiration Date	12 months from date of receipt.