

Product Datasheet

C22orf39 antibody (orb752416)

Description	C22orf39 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human UPF0545 protein C22orf39 protein (1-142aa)
Target	C22orf39
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P5X5
Dilution Range	IHC-P: 1:100-1:500
Expiration Date	12 months from date of receipt.