

Product Datasheet

C3orf38 antibody (orb752420)

Description	C3orf38 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Uncharacterized protein C3orf38 protein (80-329aa)
Target	C3orf38
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5JPI3
Dilution Range	IF: 1:100-1:500
Expiration Date	12 months from date of receipt.