

Product Datasheet

CCHCR1 antibody (orb752424)

Description	CCHCR1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Coiled-coil alpha-helical rod protein 1 protein (281-500aa)
Target	CCHCR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TD31
Dilution Range	WB: 1:1000-1:5000, IHC-P: 1:100-1:500, IF: 1:100-1:500
Expiration Date	12 months from date of receipt.