

Product Datasheet

GPRIN2 antibody (Biotin) (orb752541)

Description	GPRIN2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human G protein-regulated inducer of neurite outgrowth 2 protein (1-221aa)
Target	GPRIN2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60269
Expiration Date	12 months from date of receipt.