

Product Datasheet

KMT5C antibody (orb752600)

Description	KMT5C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Histone-lysine N-methyltransferase KMT5C protein (235-345aa)
Target	KMT5C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86Y97
Dilution Range	WB: 1:1000-1:5000, IF: 1:100-1:500
Expiration Date	12 months from date of receipt.