

Product Datasheet

LCMT1 antibody (orb752612)

Description	LCMT1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Leucine carboxyl methyltransferase 1 protein (157-357aa)
Target	LCMT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UIC8
Dilution Range	WB: 1:1000-1:5000, IF: 1:100-1:500
Expiration Date	12 months from date of receipt.