

Product Datasheet

OR13C8 antibody (orb752652)

Description	OR13C8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Olfactory receptor 13C8 protein (260-274aa)
Target	OR13C8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NGS7
Dilution Range	IF: 1:100-1:500
Expiration Date	12 months from date of receipt.