

Product Datasheet

OR9G1 antibody (Biotin) (orb752661)

Description	OR9G1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Olfactory receptor 9G1 protein (290-304aa)
Target	OR9G1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NH87
Expiration Date	12 months from date of receipt.