


## Product Datasheet

### SKA3 Antibody (orb75278)

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| <b>Description</b>         | SKA3 Antibody                                                                                |
| <b>Species/Host</b>        | Rabbit                                                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                 |
| <b>Tested Applications</b> | ELISA, WB                                                                                    |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended).<br>Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                        |
| <b>Clonality</b>           | Recombinant                                                                                  |
| <b>Uniprot ID</b>          | <b>Q8IX90</b>                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                              |



Western blot analysis of human testis ti...