

Product Datasheet

HIST1H2AG antibody (orb752828)

Description	HIST1H2AG antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Recombinant Human Histone H2A type 1 protein (27-40aa)
Target	HIST1H2AG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C0S8
Dilution Range	ICC: 1:10-1:100
Expiration Date	12 months from date of receipt.