

Product Datasheet

HIST1H2BC antibody (orb752907)

Description	HIST1H2BC antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human Histone H2B type 1-C/E/F/G/I protein (76-87aa)
Target	HIST1H2BC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62807
Dilution Range	IHC-P: 1:10-1:100
Expiration Date	12 months from date of receipt.