

Product Datasheet

Human GCIP interacting p29 protein (orb753725)

Description	Recombinant human GCIP interacting p29 protein
Reactivity	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH7.4
Storage	-20°C for 12 months as lyophilized;2-8°C for 1 month under sterile conditions after reconstitution
Tag	N-terminal His Tag
Note	For research use only
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Formula	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH _{7.4}
Purity	Greater than 95% as determined by SDS-PAGE
Source	E.Coli
MW	16.4 kDa
Uniprot ID	O95926
Expiration Date	6 months from date of receipt.