



Product Datasheet

RANBP10 Antibody (orb75387)

Description	RANBP10 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q6VN20
Expiration Date	12 months from date of receipt.



Western blot analysis of human skeletal ...



Immunohistochemical analysis of formalin...