



Product Datasheet

ENC-2 Antibody (orb75401)

Description	ENC-2 Antibody
Species/Host	Gallus
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9H0H3
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...



Immunohistochemical analysis of formalin...